# Defining Flow Rate to Reduce Initial Local-Site Adverse Reactions in Subcutaneous Immunoglobulin (SCIg) Therapy

# Background

SCIg therapy is a therapeutic alternative which an increasing number of care givers find effective. Successful self-treatment regimen requires confidence of the physician and nurse to empower independence and self-reliance for the patient. Literature and dialogue with SCIg therapy providers reveals that initial local infusion-site reactions are frequent. It is often referred to as being "expected" and "normal", but may be minimized with enhanced clinical insights which can increase patients' compliance. Immunoglobulin manufacturers recommend initial infusion volumes and flow rates which may be adjusted depending on patient outcome.

# Many patients believe that local site-reactions are to be "expected" and consider them "normal".

To minimize infusion-site reactions, knowledge of infusion system technology, technique and supplies are important factors for good tolerability and therapeutic outcome of SCIg treatments. It is also important to individualize the therapy and match it to the patient's preferences and lifestyle.

# Minimize infusion-site reactions by:



#### Having knowledge of:

- Infusion system technology
- Ancillary supplies (i.e. flow rate tubing and needle sets)



### Individualizing the therapy:

Know patient's preferences and lifestyle

#### **Contributing authors:**

Fred L. Ma, M.D., Ph.D.<sup>1</sup>, Chatarina Soderman, RN. MBA<sup>2</sup>, Andrew Sealfon<sup>3</sup>, Paul Mark Baker M.D<sup>4</sup>

<sup>1</sup>Chief Medical Officer, RMS Medical Products, <sup>2</sup>VP of Product Marketing & Clinical Affairs Associate, RMS Medical Products, <sup>3</sup> President, RMS Medical Products, <sup>4</sup>Crystal Run Healthcare

# Purpose

To assist clinicians in using a constant pressure system (CPS) to tailor SCIg therapy regimen to enhance compliance, strengthen patient confidence and control over the infusion. We present clinical strategies to understand the importance of the flow rate and its effect on outcomes.



#### **Goals for the Patient:**

- Enhance compliance
- Strengthen confidence
- Give control over infusion



# Methods

Knowing the actual flow rate (ml/hr) can be a factor in evaluating the administration and is determined by the patient simply recording the infusion time. It can be done in incremental steps from start to finish of the infusion. The flow rate can be directly calculated using the total volume divided by the infusion time.

## Using CPS, create a baseline:

# Measure 1<sup>st</sup> and 2<sup>nd</sup> half of infusion time

If both half's are equal:

- Infusion-sites are doing well
- CPS set-up is in balance with infusion-site tolerance.
- If 2<sup>nd</sup> half is longer than 1<sup>st</sup>:
- Infusion-sites are saturated
   CPS is giving feedback when saturation occurs
- Adjustments can be made to accommodate site tolerance




mode of operation (capable of providing feedback), only the CPS can be used to achieve a baseline.

Based on its unique

Flow Rate (ml) (ml/hr) = \frac{\text{Total Volume (ml)}}{\text{Infusion Time (hr)}}

# Results

Patients have confirmed actual results were very close to predicted flow rate values after several completed treatment cycles. Patients experiencing infusion difficulty reported longer infusion times.




below the 40 ml/hr flow rate line, measuring around 36ml/hr.

## Conclusion

It is important to know the actual flow rate to follow immunoglobulin manufacturers' guidelines, as well as obtaining a baseline for flow rate measurement. Clinicians can confirm the flow rate. A baseline can be created by measuring the time and volume of the infusion and compare it to the expected performance of the total infusion system including ancillary supplies. Based on the results, the clinician can make a decision to adjust the flow rate as one aspect to optimize the SClg administration.

# Peripherals can be adjusted to optimize the SCIg infusion





#### **References:**

- 1. S. Jolles, J.S Orange, A. Gardulf, M.R. Stein, R. Shapiro, M.Borte and M.Berger, 2014, Current Treatment Options with Igg for the
- Individualization of Care in Patients with PID, British Society for Immunology, Clinical and Experimental Immunology 179: 146-160.

  2. Ochoa D, Curtis C, Duff C, Riley P, Murphy E and Zampelli A. 2012. Importance of Ancillary Supplies for Subcutaneous Immunoglobulin
- Infusions: Management of the Local Infusions Site, Journal of Clinical Immunology 32 (Supplement 1): 347.

  3. Z. Panahloo, S. Jolles, 2012, Factors Affecting Infusion of High (20%) vs Lower Concentration (16%) SCIg in Primary Immunodeficiency
- Disorders, J Allergy Clin Immunol February.

  4. Sealfon A. 2012. Infusion Pump Design for SCIg: Theoretical Considerations, Journal of Clinical Immunology 32 (Supplement 1): 405.
- 5. Skoda-Smith S, Torgerson T.R och Ochs H.D. 2009. Subcutaneous Immunoglobulin Replacement Therapy in the Treatment of Patients
- with Primary Immunodeficiency Disease. Therapeutics and Clinical Risk Management 2010 (6): 1–10.

  6. R.L Wasserman, 2008, Common Infusion-Related Reactions to Subcutaneous Immunoglobulin Therapy; Managing Patient Expectation,
- Patient Preference and Adherence 2: 163-166.

  7. H. Lima and Clarke A, 2012, IgG Therapy for the Home-Based Patient: Administration and Delivery Method Considerations,
- 7. H. Lima and Clarke A, 2012, IgG Therapy for the Home-Based Patient: Administration and Delivery Method Consideration Supplement to Infusion, November/December 2012, Vol.18, No. 6.
- 8. Immunoglobulin National Society, IgNS Nursing Standard of Practice, First edition 2015.

TRS\_PO1\_DefiningFlowRate\_lgNS2016\_vA

