


# Get Started!

**FREEDOM60<sup>®</sup>**  
Syringe Infusion System



## High & Low Flow Rate Tubing Starter Kits

Available in two configurations, High for IgG and Low for antibiotics, these 20pc boxes contain a variety of our most popular flow rate tubing sets.

Starter kits are the ideal solution to minimize inventory and also provide you and your patients with the versatility to begin therapy using one flow rate, and then increasing or decreasing depending on individual needs.

**Get started today!**

### L20KT Low Flow Rate Starter Kit

Contains:

- (3) F30
- (2) F45
- (8) F60
- (4) F120
- (3) F180

### H20KT High Flow Rate Starter Kit

Contains:

- (2) F275
- (5) F600
- (5) F900
- (4) F1200
- (4) F2400

## Freedom60 Syringe Infusion System

- Advanced Technology – No batteries or electric required!
- Lowest Cost per Use – Uses standard 60cc syringes
- Reliable – No maintenance or calibration cycles required
- Tubing controls the flow rate - Wide variety of rates available
- Constant Safe Pressure – 13.5psi
- No runaway, freeflow or overdose
- Holds full pressure during and after infusion



[www.freedom60.com](http://www.freedom60.com)



845-469-2042 (USA • Chester, NY)  
info@rmsmedpro.com  
www.rmsmedpro.com