

# Finally.



***An Infusion System for Antibiotics, SCIg,  
Desferal, Chemotherapeutics, Pain  
Management and Cardiac Medications***

## **NURSING**

- Constant safe pressure—13.5psi
- No runaway, no freeflow
- Flat flow profile, less clotting
- Ambulatory, portable, lightweight
- Improves Quality of Life
- Easy to train, easy to use
- Holds pressure after infusion
- Reliable, rugged design

## **PHARMACY**

- Wide range of flow rates
- Syringes are more drug stable
- Uses standard 60cc BD® or Monoject® syringes
- Handles most specialty drugs
- Vancomycin formulations
- Desferal flow rates available

## **FINANCE**

- Reduces cost and improves healthcare quality
- Costs less than elastomerics, rentals or CADDs
- Saves on training, maintenance, batteries, nursing and drug cost (stable syringes)
- Low capital investment
- 4,000+ infusions or years of use



845-469-2042 (USA • Chester, NY)  
info@rmsmedicalproducts.com  
[www.freedom60.com](http://www.freedom60.com)

**FREEDOM60**®  
Syringe Infusion System

## Designed for the Infusion Professional.

### PERFORMANCE

- Precise infusion rates and uniform flow profiles provide consistent delivery of medication.
- Fixed flow rate sets are available from F0.5 to F2400
- Extended antibiotic admixture storage versus minibags and elastomerics.<sup>(1)</sup>

### SAFETY

- Patented luer disc syringe connector protects against set touch contamination and bolus flows.
- Highly accurate constant-force pump maintains recommended pressure at 13.5 psi.<sup>(2)</sup>
- Clear view syringe holder makes it easy to verify flow rate and volume of drug quantity.
- Protects against free flow and overdose.<sup>(3)</sup>

### CONVENIENCE

- Single or multiple syringe dosing allows extensive range of volumes (5cc to 180cc or more!).
- Easy to use, lightweight windup pump, patient just loads syringe and turns pump on.
- Completely portable, no batteries or electricity required, no cumbersome IV pole needed. Handy travel pouch available.

### LOW COST

- Lowest cost infusion system, costs less than gravity for most administrations.
- Attractively priced reusable pump (4,000+ cycle) and economical tubing sets
- Freedom60 pump is reimbursed by Medicare.
- Easy system to set-up and check. Simplifies filling, mixing and storage.
- Studies show “the syringe pump is the most cost-effective way of delivering intravenous antibiotics”<sup>(4)</sup>
- Syringes generally have greater stability and drug storage life, further reducing cost

## ORDERING INFO

### PRODUCT PART #

Freedom60 Pump	F10050
Travel Pouch	F10090

### Flow Rate Tubing Sets (50pc Box)

F0.5	F30
F1	F45
F2	F60
F3	F120
F3.8	F180
F5	F275
F8	F600
F10	F900
F15	F1200
F20	F2400

### Flow Rate Tubing Starter Kits (20pc Box)

L20KT	Low Flow Rate Tubing Starter Kit Contains: (3) F30, (2) F45, (8) F60, (4) F120, (3) F180
H20KT	High Flow Rate Tubing Starter Kit Contains: (2) F275, (5) F600, (5) F900, (4) F1200, (4) F2400

## TECHNICAL SPECS

### PUMP

#### Measurements:

Weight:	14oz (0.4kg)
Length:	12" (304mm)
Width:	4.5" (114mm)
Height:	1.6" (41mm)

#### System:

Flow rate accuracy:	+/- 8%
Operating pressure:	15psi (maximum)

### RATE CONTROLLED TUBING

#### Measurements:

Lengths:	6" to 72" (152mm - 1829mm)
----------	----------------------------

### SYRINGE

60ml BD<sup>®</sup> or Monoject<sup>®</sup> syringes only

(1) Ref Drug Stability and Compatibility Special Considerations for Home Care, H.A. Lima, R. Ph., Infusion, August 1996, pp 10-16

(2) Ref Aspects of Pressure Build-up in the Use of Electronic Infusion Devices, M.B. Stuhmeier, Anasth Intensivther, August 1997, pp 185-190

(3) Ref Sentinel Event ALERT Issue Fifteen, Joint Commission on Accreditation of Healthcare Organizations, November 30, 2000

(4) Ref Cost Comparison of Four Intermediate Intravenous Systems for Homecare Use, S. Broucek, R. Ph., American Society of Health System Pharmacists, November 1996