

FREEDOM60[®]

Syringe Infusion System 



**Antibiotics • SCIg • Desferal
Chemotherapeutics • Pain Control**



+ **Military Applications** +

- Patient Transport, KVO, TKO
- Unaffected by EMP or Electronic Warfare Signals (Hospital Ships etc...)
- Medivac or Aeromedical Use
- Lightweight, Rugged, Easy to Transport
- No IV Pole Required!
- Can be used at differing heights
- Reduce to two persons in the line of fire

- Advanced Technology – **No batteries** or electric required!
- **Lowest Cost per Use** – Uses standard 60cc syringes
- Reliable – **No maintenance** or calibration cycles required
- Tubing controls the flow rate - Wide variety of rates available
- Constant **Safe Pressure** – 13.5psi
- **No runaway, freeflow or overdose**
- **Holds full pressure** during and after infusion

MADE IN THE USA 

845-469-2042 (USA • Chester, NY)
info@rmsmedicalproducts.com
www.rmsmedicalproducts.com

**RMS**[®]
MEDICAL PRODUCTS