

Thermal Cautery Unit

Product Information

Features

Benefits

Unique Cautery "Hook"

- Controls ligation at width of 5 millimeters.
 - minimizes damage to fallopian tube.
 - increases probability of reanastomosis. (surgery to "reconnect" tubes)
- Easier to "snare" fallopian tubes.
- Ligates and cauterizes in single pass simultaneously.

Protective Teflon Cannula

- Isolates fallopian tube during transection and coagulation; shields nearby organs during procedure.
- Secures fallopian tube on "hook" when retracted.

Sterile Packaged Disposables

- Reduces risk of cross-contamination.
- Eliminates post-procedure cleaning/sterilization.

Battery Power Supply

- Low voltage for operator safety, no risk of electrical shock, burns, or electrical arcing.
- Easily transported to alternate surgical suites.
- Convenient for placement around patient area.

Other Power Supply Features:

- Line cord lockout
- Foot pedal
- Small size, lightweight

- No risk of mains current leakage.
- Convenient to operate.
- Portable, easy to carry and store.



Gyneco
24 Carpenter Road
Chester, NY 10918
800-624-9600