

## Ordering Information

### Complete Hospital Kit:

**Part #:**  
P1F1B1AFS

- RES-Q-VAC® Pump
- FSP Canister
- Soft Yankauer
- Rigid Yankauer
- 3.5"/14" 8fr Catheters
- 3.5"/14" 10fr Catheters
- 14" 14fr Catheter
- Carry Bag
- Sterile

*Also available with extra FSP Canister:*

**Part #:**  
P2F1B1AFS

- 3.5"/14" 10fr Catheters
- 14" 14fr Catheter
- Carry Bag
- Sterile

### Complete Hospital Kit Refill:

**Part #:**  
R1F1B1AFS

- FSP Canister
- Soft Yankauer
- Rigid Yankauer
- 3.5"/14" 8fr Catheters
- 3.5" 14fr Catheter
- 3.5"/14" 10fr Catheters
- 14" 14fr Catheter
- Sterile



Home Patient Kit  
Part # P1N001AFS



Nursing Home Kit  
Part # P1F1S114S



DENTAL-EVAC Kit  
Part # FSP6000

Please note: All refill kits are sold in quantities of 10.



Complete Hospital Kit\*  
Part # P1F1B1AFS  
\*shown with optional LED light.



First Responder Kit - Soft Yankauer  
Part # P1F1S110N



First Responder Kit - Rigid Yankauer  
Part # P1F1R110N

### First Responder Kit

Available with Soft or Rigid Yankauer

**Part #:**  
**Soft:** P1F1S110N  
**Rigid:** P1F1R110N

- RES-Q-VAC® Pump
- FSP Canister
- Soft Yankauer  
or Rigid Yankauer
- 3.5"/14" 10fr Catheters
- Carry Bag

### First Responder Kit Refill

Available with Soft or Rigid Yankauer

**Part #:**  
**Soft:** R1F1S110N  
**Rigid:** R1F1R110N

- FSP Canister
- Soft Yankauer  
or Rigid Yankauer
- 3.5"/14" 10fr Catheters

### Home Patient Kit:

Perfect for Ventilator or Tracheostomy Patients

**Part #:**  
P1N001AFS

- RES-Q-VAC® Pump
- Non-FSP Canister
- 3.5"/14" 8fr Catheters
- 3.5"/14" 10fr Catheters
- 3.5"/14" 14fr Catheters
- Carry Bag
- Sterile

### Home Patient Kit Refill:

**Part #:**  
R1N001AFS

- Non-FSP Canister
- 3.5"/14" 8fr Catheters
- 3.5"/14" 10fr Catheters
- 3.5"/14" 14fr Catheters
- Sterile

### Nursing Home Kit:

**Part #:**  
P1F1S114S

- RES-Q-VAC® Pump
- FSP Canister
- Soft Yankauer
- 3.5"/14" 14fr Catheters
- Carry Bag
- Sterile

### Nursing Home Kit Refill:

**Part #:**  
R1F1S114S

- FSP Canister
- Soft Yankauer
- 3.5"/14" 14fr Catheters
- Sterile

### DENTAL-EVAC™ Kit:

**Part #:**  
FSP6000

- DENTAL-EVAC™ pump
- Two canisters with FSP Filter
- Wide bore soft yankauer
- Saliva Ejector
- Carry Bag

### DENTAL-EVAC™ Kit Refill:

**Part #:**  
FSP6000R

- 2 FSP Canisters
- Wide bore soft yankauer
- Saliva Ejector

DENTAL-EVAC Pediatric Kit & Refill Also Available. Please call for more information.



# RES-Q-VAC®

Hand Held Medical Suction

## FEATURES AND BENEFITS

<b>Portable</b>	<ul style="list-style-type: none"> <li>Lightweight, 0.5 lbs</li> <li>Available whenever and wherever needed</li> <li>Compact and durable for challenging conditions</li> </ul>
<b>Manual Power</b>	<ul style="list-style-type: none"> <li>No battery maintenance or electricity required</li> <li>Works in natural disasters, or anytime the power fails</li> </ul>
<b>Reusable Pump Handle</b>	<ul style="list-style-type: none"> <li>Minimal effort required - Will not fatigue user</li> <li>Achieves high sustained air flow</li> <li>Can be ETO sterilized</li> </ul>
<b>Patented Full Stop Protection® 0.22 Micron Filter</b>	<ul style="list-style-type: none"> <li>Safeguards personnel and the general public from pathogens including: Flu, HIV/AIDS, TB, Hepatitis, SARS, and other highly infectious diseases</li> <li>Prevents fluid overflow even when the canister is filled to capacity</li> <li>Protects the pump for safe re-use</li> <li>Meets CDC, SARS and OSHA Regulations</li> </ul>
<b>Whistle-tipped Catheters</b>	<ul style="list-style-type: none"> <li>Prevents tissue damage - hospital grade</li> <li>Equipped with push connectors for ease of use</li> <li>Available sterile</li> <li>Wide range of sizes for infant through adult suctioning</li> <li>Angled rigid yankauer allows clear view of airway</li> </ul>
<b>Custom Configurations</b>	<ul style="list-style-type: none"> <li>Wide range of kits for a wide variety of markets including Hospital, EMS, nursing Homes, home care, public facilities</li> <li>Customizable for infant, pediatric and adult suctioning</li> <li>Available sterile and non-sterile</li> </ul>
<b>Powerful Vacuum</b>	<ul style="list-style-type: none"> <li>Comparable to electric wall suction - &gt;600mmHg</li> <li>Powerful enough to clear vomit, clotting blood, and mucus</li> </ul>
<b>Continuous Suctioning</b>	<ul style="list-style-type: none"> <li>Suctions until fluids or obstructions clear tubing</li> <li>Less pumping required</li> </ul>

## REPLACEMENT PARTS

<b>Catheters Non-Sterile</b>	<b>Part #:</b>	<b>Pkg:</b>
8fr. 3.5" & 14" Catheter Set	UP083514	50
10fr. 3.5" & 14" Catheter Set	UP103514	50
14fr. 3.5" & 14" Catheter Set	UP143514	50

<b>Catheters Sterile</b>		
8fr. Catheter 3.5"	835-SFR	50
8fr. Catheter 14"	814-SFR	50
10fr. Catheter 3.5"	1035-SFR	50
10fr. Catheter 14"	1014-SFR	50
14fr. Catheter 3.5"	1435-SFR	50
14fr. Catheter 14"	1414-SFR	50

<b>Yankauer Options</b>		
Soft Wide Bore 38fr 8" Sterile	R001S000S	50
Soft Wide Bore 38fr 8" Non-Sterile	R001S000N	50
Rigid Yankauer Sterile	R001R000S	50
Rigid Yankauer Non-Sterile	R001R000N	50

<b>Pump Components</b>		
Canister with FSP Filter	FSP Can	10
RES-Q-VAC® Pump - Bulk	P0000000B	10
LED Light	RQVLED	10

