

get to know

RES-Q-VAC[®]

Hand Held Medical Suction

Portable

Lightweight (0.5lbs) and compact for challenging conditions.



Manual power

No batteries or electric required – works in natural disasters, or anytime the power fails.



Reusable pump handle

Minimal effort required and can be used repeatedly.




300ml canister

Simply attaches to pump handle. Filtered and non-filtered options with gradations for easy monitoring.




Powerful vacuum

Comparable to electric wall suction at >600mmHg. Easily can clear vomit, clotting blood, and mucus.




Custom configurations

Customizable for infant through adult suctioning in sterile and non-sterile options.




Versatile

Wide range of kits for a vast variety of markets including hospital, EMS, long-term healthcare, home care, etc.




Patented Full Stop Protection[®]

0.22 Micron filter - safeguards from pathogens and infectious diseases. Prevents fluid overflow even when the canister is filled to capacity.




Whistle-tipped catheters

Hospital grade – prevents tissue damage. Equipped with push connectors for ease of use.



Continuous suctioning

Suctions until fluids or obstructions clear tubing with little pumping required.



Molded ABS plastic

Rugged, durable and operates under temperatures ranging from 0-140 degrees Fahrenheit.



Silent operation

Eliminates patient anxiety that can be triggered by complicated set-up and loud noises created by electric pumps.



Highly visible carrying pouch

Easy to locate for rapid emergency access. Provides ample storage for suctioning equipment yet is compact.

