


get to know


# HlgH-Flo™ Needle Sets

Part of the Freedom Integrated Syringe Infusion System




## Backcut needle tip design

Minimizes damage, scarring & local site reactions. Promotes faster healing.




## Flexible flying hinge

Wings have great flexibility & allow each needle to stay in place without tenting.




## Safety closure

Protects against accidental needle stick injuries.




## Custom luer

Designed to provide even flow to every site.



## Gauge options

26G needles designed for sub-q flow rates up to 40ml/hr/site & 24G for up to 300ml/hr/site.




## Length options

Needles are available in 4, 6, 9, 12, & 14mm lengths.



## Tubing length

Efficient & practical 20" length (in 26G sets) for less tangles & faster flow rates. 24" in 24G sets.



## Y-connector


Combines sets for more than 6 sites (up to 8), or to enable infusing with different length needles.




## 3M Tegaderm™

Preferred by most patients, this transparent dressing is included with every needle set.

## 26G HlgH-Flo set options & residual volume



## 24G HlgH-Flo set options & residual volume



**NOTE:** For use with a RMS Freedom Syringe Driver and RMS Freedom System accessories. Use with any other brand may cause over delivery of fluids, unknown flow rates or additional site complications such as pain, swelling and redness.

## Related Freedom Integrated System Products:



### FREEDOM60® Syringe Infusion System

For larger infusion doses, 60ml syringe



### FreedomEdge® Syringe Infusion System

For smaller infusion doses, 20ml & 30ml syringes



### RMS Precision Flow Rate Tubing™

For flow rates from F0.5 - F2400



### 24" Extension Tubing Set (FEXT)

For those who need longer tubing length



[rmsmedicalproducts.com](http://rmsmedicalproducts.com) 845-469-2042 800-624-9600

The FREEDOM60® and FreedomEdge® Syringe Infusion Systems, RMS Precision Flow Rate Tubing™, and HlgH-Flo Subcutaneous Safety Needle Sets™ are registered trademarks of RMS Medical Products and are compliant with Medical Device Directive 93/42/EEC. Tegaderm™ is a registered trademark of 3M. RMS Medical Products is ISO 13485 certified. © 2017 RMS Medical Products; All Rights Reserved. RHI\_INF\_G\_HFLbenefits\_03-10003vC